

Peripheral Venous Access Devices/Use

PAUL L BLACKBURN, RN MNA, VA-BC™

Indications for Vascular Access

- ▶ Hydration
- ▶ Medication administration
- ▶ Total or partial parenteral nutrition
- ▶ Chemotherapy
- ▶ Blood and blood component sampling/administration
- ▶ Clinical or diagnostic access
- ▶ Hemodynamic monitoring
- ▶ Arterial blood pressure monitoring/sampling
- ▶ Frequent blood sampling (DVA)
- ▶ Plasmapheresis/Aquapheresis/Hemodialysis



Device Selection



- ▶ Collaborative Consultation (multi-disciplinary team, patient/caregiver)
 - ▶ Prescribed therapy or treatment regimen
 - ▶ Anticipated duration of therapy
 - ▶ Patient assessment
 - ▶ Vascular characteristics
 - ▶ Patient's age
 - ▶ Comorbidities
 - ▶ History of infusion therapy
 - ▶ Preference for VAD location
 - ▶ Ability and resources available to care for the device

Device Selection



Peripheral Devices

- ▶ Short Peripheral Catheters (SPC)
 - ▶ Open
 - ▶ Wire assisted
 - ▶ Closed
- ▶ Midline catheters
 - ▶ Traditional
 - ▶ Multi-lumen
 - ▶ Seldinger assisted



SPC Indications

- ▶ Indications
 - ▶ Infusions projected for 6 days or less
 - ▶ Non-irritant medications and solutions
 - ▶ Non-vesicant medications and solutions



SPC Contraindications

- ▶ Contraindications
 - ▶ Thrombosis
 - ▶ Fracture/trauma
 - ▶ Infection
 - ▶ Compromised circulation
 - ▶ Mastectomy/axillary lymph dissection
 - ▶ A/V fistula or graft



Midline Catheter Indications


- ▶ Indications
 - ▶ Infusions projected for 6 days and potentially up to 4 weeks
 - ▶ Non-irritating medications and solutions
 - ▶ Non-vesicant medications and solutions



Midline Catheter Contraindications

- ▶ Contraindications
 - ▶ Fracture/trauma
 - ▶ Infection
 - ▶ Compromised circulation
 - ▶ A/V fistula
 - ▶ Chronic kidney disease patients





Decision
Process

Infusate Characteristics

Irritant

Vesicant

Osmolarity

Definitions

- ▶ **Irritant:** An agent capable of producing discomfort (burning, stinging) or pain as a result of irritation in the internal lumen of the vein with or without immediate external signs of vein inflammation.
- ▶ **Vesicant:** An agent capable of causing blistering, tissue sloughing, or necrosis when it escapes from the intended vascular pathway into the surrounding tissues.



INS List of Non-Cytotoxic Vessicants

▶ Red List

- ▶ Calcium Chloride
- ▶ Calcium Gluconate
- ▶ Contrast Media – nonionic
- ▶ Dextrose containing $\geq 12.5\%$
- ▶ Dobutamine
- ▶ Epinephrine
- ▶ Norepinephrine
- ▶ Parenteral nutrition solutions exceeding 900 mOsm/L
- ▶ Phenylephrine
- ▶ Phenytoin
- ▶ Promethazine
- ▶ Sodium bicarbonate
- ▶ Sodium chloride $\geq 3\%$
- ▶ Vasopressin

▶ Yellow List

- ▶ Acyclovir
- ▶ Amiodarone
- ▶ Arginine
- ▶ Dextrose concentrations $\geq 10\%$ to 12.5%
- ▶ Mannitol $\geq 20\%$
- ▶ Nafcillin
- ▶ Pentamidine
- ▶ Phenobarbital sodium
- ▶ Potassium $\geq \text{mEq/L}$
- ▶ Vancomycin hydrochloride

Duration of Therapy

SPC
Up to 6 days



Midline
6 days to 4
weeks



CVC
4 weeks +

Patient Assessment



Patient's age
Nutritional status



Comorbidities



History of infusion
therapy




Preference for VAD
location



Ability and
resources available
to care for VAD

SPC Site Selection


- ▶ Most likely to last duration of therapy
- ▶ Use forearm to
 - ▶ Increase dwell time
 - ▶ Decrease pain
 - ▶ Promote self-care
 - ▶ Prevent accidental removal or occlusion
- ▶ Consider veins found on the dorsal and ventral surfaces of the upper extremities
 - ▶ Metacarpal
 - ▶ Cephalic
 - ▶ Basilic
 - ▶ Median



A: Cephalic
B: Median antebrachial
C: Cephalic
D: Basilic

Midline Site Selection

- ▶ Select site in the upper arm (preferred) or
- ▶ Select site in the region of the antecubital fossa (secondary)
 - ▶ Basilic
 - ▶ Cephalic
 - ▶ Median Cubital
 - ▶ Brachial
- ▶ Avoid cannulation
 - ▶ Areas with pain on palpation
 - ▶ Open wounds
 - ▶ Infection
 - ▶ Compromised veins




Other Considerations

- ▶ Clinically Indicated
- ▶ Common modes of SPC Failure
- ▶ Bundled approach

INS Recommendations

- ▶ Remove the SPC if it is no longer included in the plan of care or has not been used for 24 hours or more.
- ▶ Remove SPC and Midline catheters when clinically indicated, based on findings from site assessment and/or clinical signs and symptoms of systemic complications.



INS Recommendations

- ▶ Potential complications
 - ▶ Pain and/or tenderness with or without palpation
 - ▶ Change in color (erythema or blanching)
 - ▶ Changes in skin temperature (hot or cold)
 - ▶ Edema
 - ▶ Induration
 - ▶ Leakage of fluid or purulent drainage from puncture site
 - ▶ Resistance when flushing
 - ▶ Absence of a blood return



Accepted but Unacceptable: Peripheral IV Catheter Failure

Peripheral IV Catheter Failure Rates

TABLE 1

Peripheral IV Catheter Failure Rate, Assorted Studies, 1990-2014

Study Type	Incidence of Failure (%)	Median	Mean
Prospective randomized controlled	36, ¹⁸ 37, ⁵⁸ 39, ⁵⁰ 40, ⁶ 45, ⁵⁹ 51, ¹⁹ 55, ²⁰ 63 ⁶⁰	43%	46%
Prospective observational	23.5, ⁶¹ 25.5, ⁶² 32, ⁶³ 36.5, ⁶⁴ 47.5, ⁶⁵ 47.6, ⁵¹ 65, ⁶⁶ 66, ⁵⁴ 69.2, ⁵⁵ 77 ⁵³	48%	49%
Retrospective	22.4, ⁶⁷ 95 ⁶⁸	58.7%	59%

5 Modes of Peripheral IV Catheter Failure

TABLE 2

The 5 Modes of Peripheral IV Catheter Failure: Prospective Randomized Controlled Studies, 1990-2014^a

Mode of Peripheral IV Catheter Failure	Range	Mean	Median
Catheter-related phlebitis	0.1%-63.3%	15.4%	9.0%
Catheter infiltration	15.7%-33.8%	23.9%	22.2%
Catheter occlusion/mechanical failure	2.5%-32.7%	18.8%	22.8%
Catheter dislodgment	3.7%-9.9%	6.9%	7.0%
Catheter-related infection	0.0%-0.44%	0.2%	0.2%

^aSummary of data from Tables 4 to 8.


Modes of IV Catheter Failure

- ▶ Overall SPC failure rate lies between 35% - 50%
- ▶ Modes of failure
 - ▶ Phlebitis
 - ▶ Infiltration
 - ▶ Dislodgement
 - ▶ Mechanical failure (eg, occlusion, leakage)
 - ▶ Site or bloodstream infection



Approach to Improve SPC Dwell Time

- ▶ Minimize phlebitis
 - ▶ Smaller-gauge catheters
 - ▶ Longer catheters
 - ▶ Catheter material
 - ▶ Stabilize catheter
 - ▶ Avoid areas of flexion
- ▶ Carefully select infusates
- ▶ Improve first attempt insertion success
- ▶ Avoid patient factors that affect vessel health



Short Peripheral Catheter Performance Following Adoption of Clinical Indication Removal

Short Peripheral Catheter Performance Following Adoption of Clinical Indication Removal

Michelle DeVries, MPH, CIC®, VA-BC ● Kathryn Strimbu

ABSTRACT

Two years following the adoption of clinical indication policies for short peripheral catheters (SPCs), a large community hospital undertook 2 extensive point prevalence reviews at 1-year intervals to study the overall outcomes associated with the SPCs. The findings were used to enhance documentation as well as staff awareness. A bundled approach was taken, focusing on insertion as well as care and maintenance needs. Consistent outcomes included at least 20% of catheters remaining functional more than 7 days and 35% more than 5 days.

Key words: clinical indication, dislodgement, dressing disruption, infiltration, peripheral IV, short peripheral catheter, vascular access

Approach to Improve SPC Dwell Time

- ▶ Catheter Infiltration
 - ▶ Avoid areas of flexion
 - ▶ Secure the catheter
 - ▶ Use stabilization
 - ▶ Carefully select catheter material
 - ▶ Limit unnecessary additional venous perforation at the time of insertion
 - ▶ Consider patient characteristics



Approach to Improve SPC Dwell Time

- ▶ Catheter occlusion and mechanical failure
 - ▶ Avoid areas of flexion
 - ▶ Stabilize the catheter
 - ▶ Optimal care and use
 - ▶ Flushing
 - ▶ Dressing care and maintenance



Approach to Improve SPC Dwell Time

- ▶ Accidental catheter removal/dislodgement
 - ▶ Stabilize the catheter
 - ▶ Use transparent dressing
 - ▶ Improved catheter use and maintenance
 - ▶ Replace dressing when loose or soiled
- ▶ Catheter-related infection
 - ▶ Don't use tape in proximity to the SPC
 - ▶ Meticulous skin prep
 - ▶ Occlusive dressing
 - ▶ Strict aseptic technique



Using a Bundled Approach to Increase SPC Dwell time

- ▶ Bundle elements
 - ▶ CHG skin prep
 - ▶ Sterile gloves
 - ▶ SPC stabilization
 - ▶ Extension set
 - ▶ Dressing materials
 - ▶ Liquid adhesive (potentially)
 - ▶ SPC insertion in forearm



References

- ▶ Resource Guide for Vascular Access. Association for Vascular Access, 2019
- ▶ Infusion Therapy Standards of Practice. Infusion Nurses Association, 2016
- ▶ Accepted but Unacceptable: Peripheral IV Catheter Failure, Journal of Infusion Nursing, Helm, R. E., et al, Volume 42, Number 3, May/June 2019
- ▶ Short Peripheral Catheter Performance Following Adoption of Clinical Indication Removal, Journal of Infusion Nursing, DeVries, M, Volume 42, number 2, Mar/Apr 2019



Thank
you